

## **Johne's Collection Container Requirements**

Starting <u>July 1, 2026</u>, we will REQUIRE fecal samples/specimens submitted for Johne's PCR to be sent in rigid plastic, leak-proof, sterile, screw cap containers. After July 1, 2026, additional processing charges will be applied to improper submissions.

<u>Please contact the NDSU VDL if you have questions or need help finding suitable container.</u> We do not promote any one brand, but examples of containers are below and can be found by suppliers such as Amazon, Consolidated Plastics, Uline, VWR, McKesson and many more.

## Acceptable:

- 25 mL or 50 mL conical tubes (free standing or U or V bottom)
- 2 oz wide-mouth jars with screw cap lid (4 oz are fine too)
- Specimen/urine collection containers (Often ~3 or 4 oz)
- Sterile is preferred, but clean and manufacturer capped container is acceptable

## Avoid:

- Condiment or dipping sauce cups as they are not sturdy enough, nor are they screw cap.
- Milk collection snap-cap tubes because gases build up and cause tubes to open and splatter.
- Sleeves because we have seen a lot of messy submissions leading to contamination between samples prior to arrival at the laboratory.
- Whirl paks, sandwich bags, grocery bags, glass containers

## Why?

- 1. Less contamination, means potentially faster and more accurate results.
- 2. Save clients money by limiting unnecessary retesting and/or culling of animals.
- 3. Ease labeling burden because they are easy to label with sample numbers (1, 2, 3 etc.) using permanent markers.
- 4. Easy to wipe down (limiting contamination) and package.
- 5. Speed up processing at the NDSU VDL which will help us turnaround results faster.
- 6. Keep testing costs down as the lab will not have to factor in as much tech time and materials.

